

# LaproClose™

## Trocar Closure Device



Intuitive, effective and safe suture placement resulting in consistent port closure, every time.

Complimentary trocar swabs permit surgeons to easily remove fluid and debris that has accumulated within the trocar, minimizing laparoscope smudging during insertion.

- **Safe**  
Reproducible, uniform closure minimizes complications and costs associated with post-operative port-site hernia.
- **Consistent**  
Provides consistent closure regardless of abdominal wall thickness
- **Versatile**  
Versatile enough to facilitate the pre-tacking of ventral hernia mesh
- **Cost Savings**  
Potential cost savings of more than \$100 when compared to other manufacturers
- **Kits**  
Suture grasper in regular or bariatric length, with suture guides available in different sizes, lengths and hole configurations to fit your needs



1625 South Tacoma Way  
Tacoma, WA 98409

Office: 888-801-9251  
Fax: 866-333-9524

customerservice@surgicalprincipals.com

Product #	Description	UOM
LCLOSE-100	LaproClose with one guide (10/12mm) and two trocar swabs (5-8mm & 10-12mm)	10/Box
LCLOSE-200	LaproClose with two guides (10/12mm & 15mm) and two trocar swabs (5-8mm & 10-12mm)	10/Box
LCLOSE-300	LaproClose XL passer with two bariatric (XL) length guides (10/12mm & 15mm) and two trocar swabs (5-8mm & 10-12mm)	10/Box

# LaproClear™

## Trocar Cleaner Swabs

Better visualization throughout the procedure.

Simple, cost-effective scope and trocar cleaning.



- **Clear Visualization**  
Quickly removes blood, fluid, and debris from scopes and trocars to maintain clear visualization throughout the procedure.
- **Improved Procedural Flow**  
Reduces lens smudging and minimizes the number of times the scope must be removed for cleaning, helping limit interruptions and maintain consistent procedural flow.
- **Gentle, Effective Materials**  
Soft microfiber cloth provides streak-free cleaning without damaging the scope lens.
- **Radiopaque, High-Absorption Design**  
Trocar swab shaft is radiopaque for visibility, with a surgical foam tip designed to absorb fluid buildup.



1625 South Tacoma Way  
Tacoma, WA 98409

Office: 888-801-9251  
Fax: 866-333-9524  
customerservice@surgicalprincipals.com

Product #	Description	UOM
LCLEAR-400	LaproClear Scope and Trocar Cleaning Kit, Disposable, Microfiber cloth and two trocar swabs (5-8mm & 10-12mm)	25/Box